



Isomass Acid Resistant Hot Plates & Accessories

ISOMASS Part Number	Part Description
A4 HOTPLATES – 110V	
ANPL-A4-1V	30x20x12cm Without Temperature Regulator
ANPL-A4-R-1V	With Rim 32x22x13cm Without Temperature Regulator
A4 HOTPLATES – 230V	
ANPL-A4-2V	30x20x12cm Without Temperature Regulator
ANPL-A4-R-2V	With Rim 32x22x13cm Without Temperature Regulator
ANPL-A4-R-2V-350°	Non-Coated High Temperature Hotplate With Rim 41x31x16cm Without Temperature Regulator (<i>Max Temp 350°C</i>)
A3 HOTPLATES – 110V	
ANPL-A3-1V-210	40x30x12cm Without Temperature Regulator 20 Amps (<i>Max Temp 210°C</i>)
ANPL-A3-1V-200	40x30x12cm Without Temperature Regulator 15 Amps (<i>Max Temp 200°C</i>)
ANPL-A3-R-1V-210	With Rim 42x32x13cm Without Temperature Regulator 20A (<i>Max Temp 210°C</i>)
ANPL-A3-R-1V-200	With Rim 42x32x13cm Without Temperature Regulator 15A (<i>Max Temp 200°C</i>)
A3 HOTPLATES – 230V	
ANPL-A3-2V	40x30x12cm Without Temperature Regulator
ANPL-A3-R-2V	With Rim 42x32x13cm Without Temperature Regulator
ANPL-A3-R-2V-350°	Non-Coated High Temperature Hotplate With Rim 53x43x17cm Without Temperature Regulator (<i>Max Temp 350°C</i>)
A2 HOTPLATES – 110V	
ANPL-A2-1V-165°	60x40x12cm Without Temperature Regulator 20A (<i>Max Temp 165°C</i>)
A2 HOTPLATES – 230V	
ANPL-A2-2V	60x40x12cm Without Temperature Regulator
ANPL-A2-R-2V-240°	With Rim 62x42x13cm Without Temperature Regulator (<i>Max Temp 240°</i>)
ANPL-A2-R-2V-350°	Non-Coated High Temperature Hotplate With Rim 75x52x17cm Without Temperature Regulator (<i>Max Temp 350°C</i>)
OTHER HOTPLATES – 110V - Without Temperature Regulator	
ANPL-30x25x4-1V	30x25x12cm
OTHER HOTPLATES – 230V - Without Temperature Regulator	
ANPL-13x13x4-2V	13x13x12cm
ANPL-15x15x4-2V	15x15x12cm
ANPL-30x25x4-2V	30x25x12cm

OVER **35** YEARS OF EXCELLENCE

ANPL-30x30x4-2V	30x30x12cm
ANPL-36x32x4-2V	36x32x12cm
ANPL-40x20x4-2V	40x20x12cm
ANPL-40x40x4-2V	40x40x12cm
ANPL-45x30x4-2V	45x30x12cm
ANPL-50x30x4-2V	50x30x12cm
ANPL-54x38x4-2V	54x38x12cm

TEMPERATURE REGULATORS FOR HOTPLATES

110V

ANREG-RNP-1V-F	Non-Programmable Temperature Regulator (P116) <i>*Available with or w/out a timer function</i>
ANREG-RP-1V-E/P1	Programmable Temperature Regulator (PXG4)

230V

ANREG-RNP-2V-F	Non-Programmable Temperature Regulator (P116) <i>*Available with or w/out a timer function</i>
ANREG-RNP-2V-A-350°	Non-Programmable Temperature Regulator (PXG4) For High Temperature Hotplate (<i>Max Temp 350°C</i>)
ANREG-RP-2V-A-350°	Programmable Temperature Regulator (PXG4) For High Temperature Hotplate (<i>Max Temp 350°C</i>)
ANREG-RP-2V-E/P1	Programmable Temperature Regulator (PXG4)

IN SITU PROBE ACCESSORIES **Compatible with 110V and 230V Temperature Regulators*

ANREG-AFF-EXT	External Display With In Situ Probe
ANPROBE-IN-SITU-1	150mm In Situ Probe
ANPROBE-IN-SITU-2	270mm In Situ Probe
ANPROBE-IN-SITU-3	450mm In Situ Probe

THERMAL PROTECTION ACCESSORIES

A4 HOTPLATE

ANTH.PROT.PTFE-A4-R-SIDE	PTFE Thermal Side Protector For A4 With Rim
ANTH.PROT.PTFE-A4-SIDE	PTFE Thermal Side Protector For A4 Without Rim
ANTH.PROT.PTFE-A4-R-BENCH	PTFE Thermal Bench Protector A4 With Rim
ANTH.PROT.PTFE-A4-BENCH	PTFE Thermal Bench Protector For A4 Without Rim

A3 HOTPLATE

ANTH.PROT.PTFE-A3-R-SIDE	PTFE Thermal Side Protector For A3 With Rim
ANTH.PROT.PTFE-A3-SIDE	PTFE Thermal Side Protector For A3 Without Rim
ANTH.PROT.PTFE-A3-R-BENCH	PTFE Thermal Bench Protector For A3 With Rim
ANTH.PROT.PTFE-A3-BENCH	PTFE Thermal Bench Protector For A3 Without Rim

A2 HOTPLATE

ANTH.PROT.PTFE-A2-R-SIDE	PTFE Thermal Side Protector For A2 With Rim
ANTH.PROT.PTFE-A2-SIDE	PTFE Thermal Side Protector For A2 Without Rim

 **Isomass Scientific Inc.**

140, 5700 – 1 Street S.W Calgary, Ab. T2H 3A9

Tel: 403-255-6631 · Fax: 403-255-6958 · Toll Free: 1-800-363-7823

www.isomass.com · sales@isomass.com

ANALAB
For analysis improvement



ANTH.PROT-PTFE-A2-R-BENCH	PTFE Thermal Bench Protector For A2 With Rim
ANTH.PROT-PTFE-A2-BENCH	PTFE Thermal Bench Protector For A2 Without Rim
STANDARD RACKS *Vessel capacity (15ml, 30ml, 60ml, 120ml) and dimensions (∅) determine the number of positions per rack	
ANRACK-A2-XXml	Standard A2 Rack (600x400x50mm)
ANRACK-A3-XXml	Standard A3 Rack (400x300x70mm)
ANRACK-A4-XXml	Standard A4 Rack (300x200x50mm)
ANRACK-A5-XXml	Standard A5 Rack (200x150x20mm)

Ordering: We are happy to accept your order(s) by email, phone or fax. We welcome all order, large or small. Visa, MasterCard and AMEX are accepted or we can open an account for your organization- please forward your enquiries to sales@isomass.com.

Quality: All Analab products are covered by a 1 year warranty for electronic and electrical elements and a 2 year warranty for graphite components (under normal use).

Can't see what you need? We offer custom plates and racks to suite your precise needs. If you do not see what you need here please contact our friendly sales team with the details of your requirements and they'll be happy to provide you with a quote, accordingly.

To request a quote, confirm availability and/or place an order please contact Isomass Scientific Inc.

Monday to Friday: 8 a.m. to 5 p.m. MTS

1-800-363-7823
sales@isomass.com
www.isomass.com